





(*) - not depicted
(**) - alternative
- Only right side depicted
- Quantities in material table for complete wheelchair

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

B75 OPTION / V500 ACTIVE				---
				Plan No.
				E00129127/A.2
Scale:	1:3	Drafted	Verified	Confirmed
Date	3.06.2026	10/06/2026	10/06/2026	
Initials	AZ	AZ	AS	





No.	Plan	Article	Qty
1	EC00336 (L)	3992230	1
	EC00337 (R)	3992231	1
2	P00050205 (L)	2670265	1
	P00048656 (R)	2670365	1
3	P00049627	26696...	2
4	P00128011 (L)	2927354	1
	P00127985 (R)	2927254	1
5	YC01489	1C47754	2
6	Clamps complete	2990321	2
7	XP03430	2997530	2
8	D4.8x19 DIN7983	1001534	2
9	M5 DIN125	1001550	2
10	M5 DIN985	1001679	2
11	M6 DIN985	1001680	2
12	X001473	1002994	2
13	Adjusting Screw	1003282	2
14	Brake Cable 1,25m	1004094	2
15	Brake lever D22	1994012	2
16	Cable Clamp Bolt M5 ---	1991919	2
		3993665**	2
17	Hand grip	3990041	2
18	Brake spring	3990133	2
19	M6 DIN985 (Black)	2994063	4
20	M6x20 ISO7380 (Black)	2994010	4
21	M6 DIN125 (Black)	2994414	4
22	Cable Hose Black Cable Hose Black	1993218	2x0,65m
		2993049**	2x0,65m
23*	Cable Hose End	1908105	4

(*) - not depicted

(**) - alternative

- Only right side depicted

- Quantities in material table for complete wheelchair

2/2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

B75 OPTION / V500 ACTIVE

Scale: 1:8	Drafted	Verified	Confirmed
Date	3.06.2026	10/06/2026	10/06/2026
Initials	AZ	AZ	AS

Plan No.

E00129127/A.2

