



Isometric view
Scale: 1:2

Material				
No.	Plan	Article	Quantity	
1	YP02247 (24CM)	2800504	1	
	YP02259 (26CM)	2800604	1	
	YP02261 (28CM)	2800704	1	
	YP02263 (30CM)	2800804	1	
	YP02249 (32CM)	2800904	1	
	YP02265 (34CM)	2801004	1	
	2	YP02159	18934	4
	3	YP02160	18935	4
4	XP00700	2990934	2	
5	XP00357	2990300	2	
6	D 6,3	1000973	4	
7	M6x12	1000700	4	
8	M6X16	1008459	4	

EP00019/A

ERGO BACK REST SAGITTA KIDS				Plan Number
				EP00019 /A
Scale: 1:2	Drafted	Verified	Confirmed	
Date	13/01/2011	13/01/2011	13/01/2011	
Initials	RM	MK	TK	