



Sheet 1/2

SEAT FRAME INOVYS 2 / NO CUSHION
 ZITFRAME INOVYS 2 / ZONDER KUSSEN

Plan No.

E000748 / G

Scale: 1:3

Drafted

Verified

Confirmed

Date

06/10/2011

19/11/2018

19/11/2018

Initials

SD

SD

PJ



E000748 / G

Material			
No.	Plan	Article	Quantity
1	X001016	1001682	2
2	X001398	1009689	1
6	Y010813	16387.. (40)	1
	Y007691	16628.. (45)	1
	Y009834	16629.. (50)	1
	YP06455	27343.. (35)	1
7	Y011640	1655100	1
8	Y011652	1655200	1
9	Y011405	16573.. (40)	1
	Y011746	16630.. (45)	1
	Y011938	16631.. (50)	1
	YP06458	27344.. (35)	1
10	Y011748	1659900	2
11	Y012123	1675800	2
14	M6 DIN125	1001564	10
15	---	1001612	2
16	M6 DIN985	1001680	1
17	M6x35 ISO7380	1006057	1
18	M6x16 DIN912	1001492	8
20	M7	1900008	1
21	M8x60 ISO7380	1900809	4
23	M8x45 ISO7380	1994473	2
24	400 N	1005602	1
28	M8x40 ISO7380	1008739	2
29	Ø32	1006111	2
30A	With velcro	1901287 (40)	1
	With velcro	1901288 (45)	1
	With velcro	1901289 (50)	1
	With velcro	2996475 (35)	1
30B	Without velcro	1913667 (40)	1
	Without velcro	1913668 (45)	1
	Without velcro	1913669 (50)	1
	Without velcro	2996476 (35)	1
31	M8 DIN125	1001556	1
32	X005575	2023600	2

Sheet 2/2

SEAT FRAME INOVYS 2 / NO CUSHION ZITFRAME INOVYS 2 / ZONDER KUSSEN				Plan No. E000748 / G
Scale: 1:3	Drafted	Verified	Confirmed	
Date	06/10/2011	19/11/2018	19/11/2018	
Initials	SD	SD	PJ	

E000748 / G