



EP00592 / B

Right side not depicted on drawing

\*Short version

Material			
No.	Plan	Article	Quantity
1	Y012555	1676204	1
	*P00029946	2647404	1
2	Y012556	1676100	1
	Y012556	1676100	1
3	YP09068 (L)	2635004	1
	YP09070(R)	2635104	1
*	P00029950(L)	2647304	1
	P00029947(R)	2647204	1
4	*XP01001	HP00194	1
	X002018	H7443	1
5	X005517	2014999	1
6	X005561	2021699	1
7	X005566 (L)	1912285	1
	X005567 (R)	1912286	1
8	XP02501	GP18200	1
9	XP03080	2996656	1
10	---	1418700	1
11	M10 DIN934	1001672	1
12	M6 DIN985	1001680	1
13	D8,4 DIN125	1003199	2
14	---	1009516	2 x 0,7 [m]
15	---	1991236	2
16	M10x19	1991432	1
17	D4x8 DIN7337A	1992805	2
18	---	2991788	0,4 [m]
	(L)450mm	2992175	1
*	(M)300mm	2992176	1
	M5x16 DIN7985	3991136	2
21	Insert M6x10	1992516	2
22	M5 DIN985	1001679	2
23	M6x20	1006999	1
24	M6x12	1007000	1
25	Sticker Vermeiren	T001002	1 <sup>2</sup>
26	M5x16	1006112	2
27	M5x12 DIN912	2994072	1
28	---	1005328	0,067 m
29	---	1009710	2 x 0,2 [m]
30	---	1009709	2 x 0,2 [m]

B66 STABILO / TRIGO				Plan No.
				EP00592 / B
Scale: 1:6	Drafted	Verified	Confirmed	
Date	02/07/2021	02/02/2022	02/02/2022	
Initials	AZ	AS	TK	