



SAFETY BELT GEMINI II 32 2020: 2996843			
No.	Plan	Article	Quantity
1	POLYPROP BAND 215/25MM PP BLACK	1009516	2,5m
2	CONDUCTOR SC /25 BLACK	1991236	6
3	BUCKLE 5POINTS GEMINI	2990693	1
4	EYELET /BUTTON FOR HEELSUPPORT 8MM	1008498	4
5	EYELET FOR HEEL SUPPORT 5 MM	1998971	2
6	POLYPROP BAND 215/40MM PP BLACK	1009517	0,35m
7	COVER SAFETY BELT GEMINI 32/40 2020	2996842	1

SAFETY BELT GEMINI II 40 2020: 2996994			
No.	Plan	Article	Quantity
1	POLYPROP BAND 215/25MM PP BLACK	1009516	3m
2	CONDUCTOR SC /25 BLACK	1991236	6
3	BUCKLE 5POINTS GEMINI	2990693	1
4	EYELET /BUTTON FOR HEELSUPPORT 8MM	1008498	4
5	EYELET FOR HEEL SUPPORT 5 MM	1998971	2
6	POLYPROP BAND 215/40MM PP BLACK	1009517	0,5m
7	COVER SAFETY BELT GEMINI 32/40 2020	2996842	1

EP00620 / A

U52 SAFETY BELT GEMINI 2020				Plan No.
				EP00620 / A
Scale: 1:10	Drafted	Verified	Confirmed	
Date	21/09/2021	21/09/2021	21/09/2021	
Initials	OY	AZ	TK	