



Drawing shows only right side

EP00378 / D

ADJUSTABLE BACKREST WITH FOLDABLE HANDLE TRIGO				Plan No. EP00378/D
Page 1/2				
Scale: 1:5	Drafted	Verified	Confirmed	
Date	14/02/2019	08/10/2021	08/10/2021	
Initials	AS	MK	TK	

Material (31-36)			
No.	Plan	Article	Quantity
1	YP06265	2785265	1
2	YP07062	2744565	1
2**	YP07060	2744665	1

Material (35-40)			
No.	Plan	Article	Quantity
1	YP06265	2785265	1
2	YP07060	2744665	1
2**	YP07061	2744765	1

Material (40-45)			
No.	Plan	Article	Quantity
1	P00025189	2643065	1
2	YP07060	2744665	1
2**	YP07061	2744765	1

Material (45-50)			
No.	Plan	Article	Quantity
1	P00025189	2643065	1
2	YP07061	2744765	1
2**	YP07063	2744865	1

Material			
No.	Plan	Article	Quantity
3	M6x30	1992721	1
4	Washer M6	2994068	1
5	M6	2994059	3
6	XP01560	2993561	4
7	M6x40	2994066	2
8	Fold. push handle	EP00246	1
9	Y004353	F000411	1
10	XP02043	2993671	1
11	M4	1001678	1
12	M4x16	2994070	1
13	M6x10	2994127	3
14	M6X10	2994284	2
* 15	Cap D22	1990648	1

* When unused folding handle / Gdy nie uzywana raczka skladana.

** Gdy nie stosujemy poz. 8

**If no. 8 isn't used.

EP00378 / D

ADJUSTABLE BACKREST WITH FOLDABLE HANDLE TRIGO				Plan No.
				EP00378/ D
Scale: 1:5	Drafted	Verified	Confirmed	
Date	14/02/2019	08/10/2021	08/10/2021	
Initials	AS	MK	TK	